

STERILISING FLUID



Legend Sterilising Fluid easily and effectively disinfects hard surfaces, offering protection against bacteria, yeast, fungi and viruses. Suitable for use on hard surfaces, floors, walls and utensils. Concentrated solution, dilute with clean water. **Passes EN1276, EN13697, EN14476.**

Pack size: 2x5L Code: LEGM-AS0140-2N



How To Use

Clean surfaces with Legend Hard 'N' Fast or Legend BC40 then rinse with clean water. Dilute before use - Dilute at 3% i.e. 30ml concentrate plus 970ml clean water, per litre.

Apply the pre-diluted Sterilising Fluid through a trigger spray, soak application or by mopping.

Note: For full virucidal and bacterial kill, a contact time of 15 minutes should be allowed. Rinse with clean water ensuring any residue is removed.



DO NOT USE WHERE ACID BASED CHEMICALS ARE PRESENT AS THIS MAY PRODUCE A TOXIC GAS

HEALTH & SAFETY

Please refer to the MSDS for full health and safety information.



DO NOT INGEST
DO NOT MIX WITH OTHER CHEMICALS
KEEP OUT OF REACH OF CHILDREN